WoundLyser 2M

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INTENDED PURPOSE

The WoundLyser 2M is a medical water filter with a hollow-fibre membrane. It provides sterile-filtered tap water for irrigating chronic wounds and is screwed onto a shower hose with a G½" thread in place of the hand shower. The hygienic design of the filter and the coated spray disc help to minimise the risk of microbial growth. Wounds should only be irrigated with sterile-filtered water after consulting the treating doctor, nursing service or wound specialist.

IMPORTANT INFORMATION BEFORE USE

Protect the filter: Do not drop it, do not immerse it in wash water, do not open or damage it. Only use undamaged filters.

- **Read these instructions:** Read these instructions carefully before installation and use.
- **Medical consultation:** Only irrigate wounds after consulting the treating doctor, nursing service or wound specialist.
- **Hygiene:** To prevent microbial contamination from splashes or contact with wound material, carry out hygiene steps B1–B4 before and after every wound cleansing (p. 4).
- **Check the temperature:** Before use, check the water temperature on healthy skin, e.g. on your hand.
- **Low pressure:** Regulate the irrigation pressure at the mixer tap and start with a very low setting.
- **No low-pressure system:** Do not connect to low-pressure systems such as unpressurised boilers. Max. 5 bar, max. 60 °C.
- **Hoses:** Only use shower hoses in perfect condition; replace porous or defective hoses.
- **Service life:** Max. 2 months (62 days) from first use, then dispose of it. Do not use if the packaging is damaged or the use-by date has passed.

Storing the filter in its original packaging: keep dry and protected from moisture at 10–50 °C. Use within 24 months of the production date (barcode MMYXXXXX). Do not use if the packaging is damaged.

→ p. 2

INSTALLATION AND USE → p. 3

- A1** Unscrew the existing hand shower from the shower hose. Only use hoses in perfect condition; replace porous or defective hoses.
- A2** Screw the WoundLyser hand-tight onto the shower hose (G½").
- A3** Turn on the water supply and check for leaks: no water may escape from the coupling thread.
- A4** Carry out hygiene steps B1–B4. Check the water temperature on healthy skin (e.g. your hand), set a very low irrigation pressure at the mixer tap to begin with and irrigate the wound.
- A5** Do not immerse the filter in wash water or drop it onto hard surfaces. Use the shower head holder to hang it securely.
- A6** Never open or damage the filter. Only use undamaged filters.

HYGIENE, SERVICE LIFE, DISPOSAL → p. 4 · B1–B4 before and after every wound cleansing

- B1** Let water run through the filter for approx. 10 seconds.
- B2** Then dry the spray disc and the outside of the filter housing with a clean, dry cloth.
- B3** Next, disinfect the spray disc and housing with a commercially available alcohol-based disinfectant wipe; observe the contact time of the disinfectant.
- B4** Then let water run through the filter for approx. 10 seconds.
- B5** The service life is max. 2 months (62 days) after first use. Enter the date in the log (p. 20).
- B6** At the end of its service life, dispose of the used WoundLyser with household waste.

Order a new WoundLyser

WoundLyser 2M (item no. 1384) – scan the QR code or order online:
sanquell.de/produkt/wunddusche-wundreinigung-ambulant-log7

Video

How to use – scan the QR code

TECHNICAL DATA

Filter type	medical water filter with hollow-fibre membrane
Pore size	0,2 µm
Retention	7 log, <i>Brevundimonas diminuta</i>
Service life	max. 2 months (max. 62 days) ¹
Flow rate	10 l/min (5 bar, 30 °C)
Max. operating pressure	5 bar ³
Max. operating temperature	60 °C (70 °C ≤ 30 min over the service life)
Chlorine resistance	≤ 10 ppm ²
Connection	G½"
Dimensions	L 220 mm · Ø 44 mm

¹ Depending on water quality.

² Continuous dosing of ≤ 10 ppm over the service life; short-term (1 h) high dose (400,000 ppm) for chemical disinfection.

³ The WoundLyser must not be connected to low-pressure systems (e.g. boilers).

FREQUENTLY ASKED QUESTIONS (FAQ)

- ? **What is the WoundLyser for?**
For irrigating chronic wounds with sterile-filtered tap water – after consulting the treating doctor, nursing service or wound specialist.
- ? **Will the filter fit my shower?**
Yes, if the shower hose has a standard G½" thread. Not suitable for low-pressure systems such as unpressurised boilers.
- ? **How long can I use one filter?**
Max. 2 months (62 days) from first use, depending on water quality. Record the date in the log on p. 20.
- ? **Why flush and disinfect before and after every use?**
Splashes or contact with wound material can contaminate the outside of the filter. Steps B1–B4 prevent this.
- ? **Is the filter sterile?**
The filter is supplied non-sterile. It is the water flowing through the membrane that is sterile-filtered – which is why external hygiene is so important.
- ? **Which temperature and pressure?**
Always check the temperature on healthy skin, max. 60 °C. Regulate the pressure at the mixer tap and start very low.
- ? **I dropped the filter – what now?**
Only use undamaged filters. If there are cracks or visible damage, stop using the filter and replace it.
- ? **Can I clean or reuse the filter?**
No. Do not open or reprocess the filter. Dispose of it with household waste after its service life.

MANUFACTURER AND DISTRIBUTOR (GPSR)

Manufacturer: Aqua free GmbH
Winsbergring 31 · 22525 Hamburg · Germany
info@aqua-free.com · www.aqua-free.com

Distributor: Sanquell GmbH
Sieker Berg 9 · 22962 Siek
hallo@sanquell.de · www.sanquell.de

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Fabricante: Aqua free GmbH
Distribuidor: Sanquell GmbH